

# Marquis™

## Fresh Scented Granular Presoak Concentrate

Marquis™ is a citrus scented, granular presoak that comes in a unique canister for safe and convenient use with the Freedom System™ dispenser. Its alkaline/surfactant blend penetrates soils quickly to remove tough food soils from silver and stainless flatware.

Specifications	
Color	Blue
Scent/Odor	Citrus
Character	Granular powder
Foam	Low
pH Use Dilution	9.5 - 11.0
Phosphorus Content	3.4%
Approvals	cRC

  

Pack Sizes	
5898930	4/7 pound case



### Directions For Use

**Equipment:** Recommended dispensing system.

1. Dispense product solution into bus tub. Consult your sanitation specialist for proper use dilution.
2. Place stainless ware and/or silverware in solution for several minutes.
3. Remove, rinse and wash in dishmachine in normal fashion.

**CANISTER SET UP:**

1. Remove the small outer cap with flip tab from canister. Do not remove larger inner cap.
2. Place the canister into the Freedom System™ dispenser, open end down. Use handle when transporting canister.
3. DO NOT remove concentrate canister until completely empty.

4. To remove empty canister, lift it from the dispenser carefully by the handle and turn upright to avoid spills. Dispose of canister immediately. DO NOT reuse empty canister.
5. Place new concentrate canister into dispenser.

**NOTE:** Dispenser WILL NOT activate unless a canister is in unit.

### Product Demo



Scan the QR code with your mobile device to view an instructional video featuring this product!

### Safety Reminder

Read and understand the label and Safety Data Sheet (SDS). These provide a description of hazards, precautions, personal protective equipment, first aid, 24 Hour Emergency lines, storage data and other safe handling information. The SDS is available on our website at [www.uschemical.com](http://www.uschemical.com) and is at your fingertips by scanning the SDS QR Code on the label.

